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## Referências

1. Achilli A, Tura F, Euwe E. Immediate/early function with tapered implants supporting maxillary and mandibular posterior fixed partial dentures: Preliminary results of a prospective multicenter study. *J Prosthet Dent* 2007; 97:S52-S58.
2. Attard NJ, David LA, Zarb GA. Immediate loading of implants with mandibular overdentures: One year clinical results of a prospective study. *Int J Prosthodont* 2005;18:463-470.
3. Cehreli MC, Uysal S, Akca K. Marginal bone level changes and prosthetic maintenance of mandibular overdentures supported by 2 implants: A 5-year randomized clinical trial. *Clin Implant Dent Relat Res*. 2010 Jun 1;12(2):114-21. Epub 2009 Feb 13.
4. Eliasson A, Blomqvist F, Wenneberg A, Johansson A. A retrospective analysis of early and delayed loading of full-arch mandibular prostheses using three different implant systems: Clinical results with up to 5 years of loading. *Clin Impl Dent Relat Res* 2009;2: 134-148.
5. Fischer K, Bäckström M, Sennerby L. Immediate and early loading of oxidized tapered implants in the partially edentulous maxilla: A 1-year prospective clinical, radiographic, and resonance frequency analysis study. *Clin Impl Dent Relat Res* 2009;2:69-80.
6. Friberg B, Jemt T. Rehabilitation of edentulous mandibles by means of five TiUnite implants after one-stage surgery: A 1-year retrospective study of 90 patients. *Clinical Implant Dentistry and Related Research* 2008;1:47-54.
7. Olsson M, Urde G, Andersen JB, Sennerby L. Early loading of maxillary fixed cross-arch dental prostheses supported by six or eight oxidized titanium implants: Results after 1 year of loading, case series. *Clin Impl Dent Relat Res* 2003;1(suppl):81-87.
8. Turkyilmaz I, Tözüm TF, Turner C, Ozbek EN. A 2-year clinical report of patients treated with two loading protocols for mandibular overdentures: Early versus conventional loading. *J Periodontol* 2006;77:1998-2004.
9. Turkyilmaz I, Avci M, Huran S, Ozbek EN. A 4-year prospective clinical and radiological study of maxillary dental implants supporting single-tooth crowns using early and delayed loading protocols. *Clin Impl Dent Relat Res* 2007;4:222-227.
10. Vanden Bogaerde L, Rangert B, Wendelhag I. Immediate/Early function of Brånemark System® TiUnite implants in fresh extraction sockets in maxillae and posterior mandibles: An 18-month prospective clinical study. *Clin Impl Dent Relat Res* 2007;1(suppl): S121-S130.
11. Vanden Bogaerde L, Pedretti G, Dellacasa P, Mozzati M, Rangert B, Wendelhag I. Early function of splinted implants in maxillas and posterior mandibles, using Brånemark System TiUnite implants: An 18-month prospective clinical multicenter study. *Clin Impl Dent Relat Res* 2007;1(suppl): S121-S130.
12. Villa R, Rangert B. Immediate and early function of implants placed in extraction sockets of maxillary infected teeth: A pilot study. *J Prosthet Dent* 2007;97: S96-S108.
13. Zembic A, Glauser R, Khraisat A, Hämmerle CH. Immediate vs. early loading of dental implants: 3-year results of a randomized controlled clinical trial. *Clin Oral Implants Res*. 2010 May;21(5):481-9.